



Overview of Physigraph's Web Application

The most evolutionary program
on the market with more than
3500 image of exercises and more than
1500 animations and videos.



From the starting page
insert your user name and
password.



Insert client First and Last name
with e-mail address to be sent
the fitness program with one click.

Manage Protocols - Windows Internet Explorer

http://www.physigraph.net/index.php?do=exercices

PHYSIGRAPHE Clients Exercises Protocols Archive Administration

Logged as Sylvain Lemaire [4] · Preferences · Help · Logout

SEARCH BROWSE PROGRAM

FIND EXERCISES

Find an Exercise Select a keyword Standard Protocols: My PTCL All

Cardio
 Flexibility
 Plyometry
 Strengthening
 Warm-up

Muscles

Heart

From the "Find Exercises" page choose the activity where you want to get the exercises.

[Clients](#) · [Exercises](#) · [Protocols](#) · [Archive](#) · [Preferences](#) · [Help](#) · [Logout](#)

Copyright © 1996-2007 All Rights Reserved · [Terms of Use](#) · [License Agreement](#)

PHYSIGRAPHE Clients Exercises Protocols Archive Administration

Logged as Sylvain Lemaire [1728] · Preferences · Help · Logout

SEARCH BROWSE PROGRAM

FIND EXERCISES

Find an Exercise Select a keyword Standard Protocols: My PTCL All

Cardio	Muscles	Equipments	Positions	Others
<input type="checkbox"/>	<input type="checkbox"/> Abdominals	<input type="checkbox"/> Ball	<input type="checkbox"/> Declined	<input type="checkbox"/> Aquatic
<input type="checkbox"/>	<input type="checkbox"/> Abdominals oblique	<input type="checkbox"/> Bar	<input type="checkbox"/> Hang	<input type="checkbox"/> Weight lifting
<input type="checkbox"/>	<input type="checkbox"/> Abductors	<input type="checkbox"/> Bench	<input type="checkbox"/> Inclined	<input type="checkbox"/> Yoga
<input type="checkbox"/>	<input type="checkbox"/> Adductors	<input type="checkbox"/> Bosu	<input type="checkbox"/> Lying	
<input type="checkbox"/>	<input type="checkbox"/> Back	<input type="checkbox"/> Dumbbell	<input type="checkbox"/> Neutral Grip	
<input type="checkbox"/>	<input type="checkbox"/> Biceps	<input type="checkbox"/> Elastic	<input type="checkbox"/> On 1 foot	
<input type="checkbox"/>	<input type="checkbox"/> Buttock	<input type="checkbox"/> Free	<input type="checkbox"/> On 2 feet	
<input type="checkbox"/>	<input type="checkbox"/> Calf	<input type="checkbox"/> Heavy ball	<input type="checkbox"/> On 2 knees	
<input type="checkbox"/>	<input type="checkbox"/> Forearm	<input type="checkbox"/> Kettlebell	<input type="checkbox"/> On all fours	
<input type="checkbox"/>	<input type="checkbox"/> Hamstrings	<input type="checkbox"/> Machine	<input type="checkbox"/> On one knee	
<input type="checkbox"/>	<input type="checkbox"/> Hand	<input type="checkbox"/> Misc	<input type="checkbox"/> On the back	
<input type="checkbox"/>	<input type="checkbox"/> Hip	<input type="checkbox"/> Partner	<input type="checkbox"/> On the side	
<input type="checkbox"/>	<input type="checkbox"/> Legs	<input type="checkbox"/> Plate	<input type="checkbox"/> On the stomach	
<input type="checkbox"/>	<input type="checkbox"/> Neck	<input type="checkbox"/> Pulley	<input type="checkbox"/> Pronation Grip	
<input type="checkbox"/>	<input type="checkbox"/> Pectorals	<input type="checkbox"/> Roller	<input type="checkbox"/> Seated	
<input type="checkbox"/>	<input type="checkbox"/> Quadriceps	<input type="checkbox"/> Roller board	<input type="checkbox"/> Stabilization	
<input type="checkbox"/>	<input type="checkbox"/> Quadriceps - buttock	<input type="checkbox"/> Sliding Disc	<input type="checkbox"/> Stand	
<input type="checkbox"/>	<input type="checkbox"/> Rotators	<input type="checkbox"/> Stairs	<input type="checkbox"/> Supination Grip	
<input type="checkbox"/>	<input type="checkbox"/> Shoulder	<input type="checkbox"/> Step		
<input type="checkbox"/>	<input type="checkbox"/> Trapezius	<input type="checkbox"/> Swiss Ball		
<input type="checkbox"/>	<input type="checkbox"/> Triceps	<input type="checkbox"/> Towell		
<input type="checkbox"/>	<input type="checkbox"/> Trunk	<input type="checkbox"/> TRX		
		<input type="checkbox"/> Wooble Board		

From The Physigraph search engine search by the muscle groups, equipments, positions and others.

[Clients](#) · [Exercises](#) · [Protocols](#) · [Archive](#) · [Preferences](#) · [Help](#) · [Logout](#)

Copyright © 1996-2008 All Rights Reserved · [Terms of Use](#) · [License Agreement](#)

BROWSE AN EXERCISE

The Search Yielded 203 Exercise(s) Use the checkboxes to select exercises and the eye to view details. Select: All / None

Assign to Program Search Again



Click on the left top corner to select the exercise you want and assign it to the program.

To select other exercises go back to the search page.

BROWSE A PROGRAM

Print Size: Small Large
 Language: English French Spanish
 Expiration of online report:
 Show in B&W: Or Color

Note(s)

Do the program every 2 days

Update the note

When you've chosen your exercises, you can place them in order and then print them in black and white.

Save this Program Print this Program Done with this Program

Created by: · Date: 2012-02-2 [0.00 Min.]



To complete the program click on the eye on the top left corner of the image.

Flip (Mirror Images)

Start

End

Instructions

My Instructions
 Physigraph
 Program's instructions

- Stand with your feet wide apart, keep your back straight and abs tight - Go down on your right and on your left side

Set as my instructions

TOTAL

Volume:	<input type="text" value="0"/> KG	Temps:	<input type="text" value="18.00"/> Min.
Tension:	<input type="text" value="0.00"/> Min.	Rep:	<input type="text" value="180"/>

The information here is automatically saved when the program is saved

Alternate:

1RM: KG

Fr/Week:
 Sets:
 Rep.:
 Tpo:
 Tempo:
 Rest:
 %:
 Load:

From this page you can write all the information needed to perform the exercise and the program: like the weekly frequency, number of sets, number of repetitions, the tempo and the rest. The time of the workout, number of repetition, the time under tension and the volume will be calculated automatically.

You can add new sets to change the quality of the workout. The load could also be calculated from the 1 RM (1 repetition max).

The Physigraph's Web Application offers alot more, like animations, the creation of protocols, a complete management of your Fitness Trainers for administrators of fitness centers.

With this bank of over 3500 exercise images and over 1500 animations and videos, Physigraph's Web Application wants to be the most evolutionary program on the market to give to Fitness Trainers the best fitness tool ever offered.